



Instructions for Use



English

Corowell COVID-19 Symptom Screening Rapid Test

User Manual v1.1

MKG Consulting UG
Hauptstrasse 96
53474 Bad Neuenahr-Ahrweiler
Germany

Telephone: +49 2641 908 28 74
E-Mail: info@corowell.com

1 Intended Use

The *Corowell COVID-19 Symptom Screening Test* is intended to objectively testing of the sense of smell. It enables the identification of subjects whose olfactory perception is deviating, restricted or newly got lost.

The *Corowell COVID-19 Symptom Screening Test* is a device for home use.

The *Corowell COVID-19 Symptom Screening Test* can be performed self-administered by the subject.

The target population is adults over 18 years of age.

2 Contra Indications

The Corowell COVID-19 Symptom Screening Test is contraindicated for subjects with hereditary, pre-existing or long-term impaired olfactory perception.

3 Instructions

Please read this manual before using the Corowell COVID-19 Symptom Screening Test for the first time.

3.1 Warnings

- The Corowell COVID-19 Symptom Screening Test has been developed and tested for use by a single user. Do not use the Corowell COVID-19 Symptom Screening Test with other people, but only alone. However, you can use the same app to test multiple users in a row.
- Please do not store Corowell Tickets near water or in damp places.

3.2 Precautions

- The Corowell COVID-19 Symptom Screening Test should only be used in an environment where there are few strong odors so that the user can easily and not distractedly detect the odors on the Corowell Ticket.
- The Corowell COVID-19 Symptom Screening Test should only be used, if the user has a working Internet connection.

4 General Functions

The Corowell COVID-19 Symptom Screening Test System is designed to support the user during the process so that the user can successfully test at home or in the field.

The Corowell COVID-19 Symptom Screening Test System consists of

- Corowell Application (iOS and Android)
- Corowell Ticket



Figure 1 - Corowell COVID-19 Symptom Screening Test System Components

5 Using the Corowell System

The detailed use of the Corowell system is described in the following chapters.

5.1 Loading the App Software for the First Time

To download the app software to the user's smartphone for the first time, the user can

Search for Corowell in the respective App Store / Google Store and download and install the app from there

On the back of the Corowell ticket, scan the QR Code as shown in Figure 2. From there, the user is directed to the App Store and can download and install the App software without having to search.

Scan this QR Code to download the App



Figure 2

5.2 Start Screens

When you open the test for the first time, the following screens appear, which serve as an introduction to the screening test. They only appear when the app is opened for the first time and not later (Figure 3).




Figure 3

5.3 Ticket Scanning Screen

After the Start Screens the "Ticket Scanning Screen" appears. First you have to allow the app to access your camera (Figure 4).

5.4 Select Language

When you start the Corowell app for the first time, the software automatically adopts the language selection screen set in your smartphone. You can change this at any time in the settings (Figure 4).




Figure 4

5.5 Ticket Scanning

To scan a new ticket each time, you must first carefully open the glued ticket and then scan the QR code on the inside using the app (Figure 5). To do so, simply press the "+" sign in the blue crest at the bottom of the ticket scanning screen (figure 6).



Figure 5




Figure 6

5.6 Select Scent Screen

After you have scanned the QR code on the inside of the ticket, the Scent Selection Screen will open (Figure 7). There you have to select the scent you have detected on the green scent area (Figure 5) "after" scratching it with your fingernail or a coin.

Then please click "Confirm" at the bottom of the screen (Figure 7).




Figure 7

5.7 COVID-19 Questionnaire Screen

Here the seven (7) questions must answer ALL truthfully (Figure 8). Only if you have answered YES or NO to ALL questions and the slider has been adjusted accordingly, you can continue by clicking Confirm.




Figure 8

5.8 Enter Name and Data Privacy

Next, the screen opens as shown in Figure 9, where you must enter your name (or the name of the person who actually took the test, e.g. your children).

This is very important because later, when you present your Corowell PASS for third-party verification, you may be asked for your ID to confirm your identity.

You must also accept the privacy policy, which you can read in detail at any time either via the link (Figure 9) or in the settings.




Figure 9

5.9 Corowell PASS Screen

If you have "passed" the Corowell COVID-19 Symptom Screening Test, the Corowell PASS Screen (Figure 10) will open. Here you will see the name you entered earlier and the date and time you took the test. This is needed, so that when your Corowell PASS is being tested by a third party, the time and date at which you passed is important.


Operators of hotels, restaurants, airlines, etc. will only accept Corowell passes that are not older than 6, 12, 24, 48 or 72 hours. These acceptance criteria will be communicated to you on site.

The sudden loss of smell is one of the key early symptoms of COVID-19, but it is not a universal symptom. The ability to identify the scents provided is not a guaranty of a negative COVID-19 status.



Figure 10

5.10 Corowell FAILED Screen



If you have not passed the Corowell test, this screen will appear. If you have the impression that you "only" did something wrong, you can always repeat the test with a new ticket. If you really can't smell, or if you answered YES to one of the first 5 questions or NO to the last question, we would recommend that you have yourself examined by a healthcare professional for a possible COVID-19 infection. It is recommended to follow up by an additional COVID-19 test method, which is authorized by the local government.

6 Error Message „Sorry, we are unable to process your Ticket“

This screen (Figure 11) appears, if you try to use your ticket a second time or if you have a copy or a fake of the original Corowell ticket.




Figure 11

7 Settings

In the settings (Figure 12), you can open this user manual and read the FAQs.

Here you will also find the privacy policy (in your respective language) and can visit our corowell.com website.

In the field "About" the official medical device will find information about our product (Figure 13).




Figure 12


Figure 13

8 Service and Warranty

We offer you a comprehensive global warranty protection that comes into effect on the date of purchase. Please ask your seller for details and warranty period. This warranty covers all defects in materials and/or workmanship.

The warranty does not cover damage resulting from improper use or care, exposure to chemicals, immersion in water or inappropriate exposure. Damage caused by third parties or unauthorized service centers will void the warranty.

9 Symbols



Various symbols are used in the above information. These are among explained below:

	Medical Device
	Reference Number
	Serial Number
	Manufacturing Date
	Manufacturer
	Keep dry
	Single Use Device
	Protect from Sunlight
	Temperature Limits
	See User Manual
	European Conformity marking